

- APPROPRIATE SETTINGS FOR THE INSERTS
- SETTAGGI APPROPRIATI PER GLI INSERTI
- KORREKTE EINSTELLUNGEN FÜR DIE INSTRUMENTE
- RÉGLAGES APPROPRIÉS POUR LES INSERTS
- AJUSTES APROPIADOS PARA LOS INSERTOS
- KORREKTA INSTÄLLNINGAR FÖR SPETSARNA

→ MULTIPIEZO PRO → MULTIPIEZO touch → MULTIPIEZO white



function	endo	perio/scaler	restorative/prosthetics
<b>SCALING</b>			
S1-5 S1, S2, S3, S4, S5, S6, S7, S8		1-2* 4-6	
<b>PERIO</b>			
P1, P2, P3, P4, P10, P11, P12, P13, P14		1-3	
<b>IMPLANT CLEANING</b>			
ICS		1	
<b>ENDO</b>			
E1, E2	1-6		
<b>ENDO REVISION</b>			
ER1 ER2, ER3, ER4, ER5		1-6 1-3	
<b>ENDO RETRO</b>			
R1, R2, R3, R4, R5, RD3, RD4, RD5		1-3	
<b>CROWN PREP</b>			
DB1			3-6
<b>CAVITY MARGIN</b>			
CM1, CM2, CM3, CM4			1-5
<b>RESTORATIVE</b>			
D1, D2, D4 D3 ME1, ME2, ME3**			1-6 1-5 6

\* it is recommended the use of „soft mode“ function

\*\* only available in the United States

- APPROPRIATE SETTINGS FOR THE INSERTS
- SETTAGGI APPROPRIATI PER GLI INSERTI
- KORREKTE EINSTELLUNGEN FÜR DIE INSTRUMENTE
- RÉGLAGES APPROPRIÉS POUR LES INSERTS
- AJUSTES APROPIADOS PARA LOS INSERTOS
- KORREKTA INSTÄLLNINGAR FÖR SPETSARNA



function	endo	perio/scaler	restorative	power	power
<b>SCALING</b>					
S1-5 S1, S2, S3, S4, S5, S6, S7, S8		1-2 * 4-5		1-2 4-6	1-2 3-4
<b>PERIO</b>					
P1, P2, P3, P4, P10, P11, P12, P13, P14		3		1-3	1-2
<b>IMPLANT CLEANING</b>					
ICS		1		1-3	1-2
<b>ENDO</b>					
E1, E2	1-5			1-2	1
<b>ENDO REVISION</b>					
ER1 ER2, ER3, ER4, ER5		1-5 1-3		3-6 1-3	2-4 1-2
<b>ENDO RETRO</b>					
R1, R2, R3, R4, R5, RD3, RD4, RD5		1-3		1-4	1-3
<b>CROWN PREP</b>					
DB1			3-5	1-3	1-2
<b>CAVITY MARGIN</b>					
CM1, CM2, CM3, CM4			1-4	1-5	1-3
<b>RESTORATIVE</b>					
D1, D2, D4 D3 ME1, ME2, ME3 **			1-5 1-4 5	4-6 4-6	3-4 3-4

\* it is recommended the use of „soft mode“ function  
\*\* only available in the United States